


Multifocal contact lens fitting guides.



CooperVision®



With over 90 times more prescription options than the other main manufacturers combined, CooperVision® has by far the largest multifocal contact lens range.¹

We understand the importance of ease and accuracy during the fitting process.

These fitting guides have everything you need to help make it easy to achieve multifocal fitting success.

CooperVision® advises **use of the +1.00 blur test** for assessing sensory dominance and to achieve a successful fitting of soft multifocal lenses.

1. CooperVision Data on File 2023. Based on prescription option combinations (sph, cyl, axis & add) available across all soft lenses in multifocal, SIH and hydrogel from J&J, Alcon, B+L and CVI in UK, France, Germany and Italy Feb 2023. Includes stocked and made to order lenses in daily, 2 weekly and monthly disposable options. Cosmetic & Photochromatic CIs not included. Multiple base curve variants not included.



MyDay® multifocal fitting guide

Featuring the CooperVision® Binocular Progressive System™ and Aquaform® Technology

- Simple to fit and an easy progression as the near requirement increases.¹
- 98% of patients successfully fitted with two pairs of lenses or fewer.^{1,2}
- 83% of patients successfully fitted with the first pair of lenses.³
- Aquaform® Technology for incredible all-day comfort.



Initial lens selection

- Step 1** – Using up-to-date spectacle prescription, determine ocular dominance using the **+1.00D blur** method.
- Step 2** – Determine spherical equivalent distance power (corrected for vertex distance).
- Step 3** – Select distance sphere power for each eye with ADD powers as indicated below.

FITTING GUIDE		
Spectacle Rx	Dominant Eye	Non-Dominant Eye
ADD	DESIGN	DESIGN
+0.75 to +1.25D	LOW	LOW
+1.50 to +1.75D	LOW	MED
+2.00 to +2.50D	LOW	HIGH

Vision assessment

- For best results, allow wearer to experience vision outside the testing room for 10-15 minutes.
- Check vision with both eyes open and room lights on.
 - For **distance vision**, assess in surrounding environment under normal lighting conditions.
 - For **near vision**, assess using a mobile phone or other reading material.
- If vision acceptable, dispense trial lenses.
- If vision not acceptable, follow the lens optimisation steps described to the right.

Lens optimisation



Have patient keep both eyes open and optimise using handheld lenses or a flipper. **DO NOT USE A PHOROPTER.**

DO NOT CHANGE ADD POWER.

DISTANCE VISION ENHANCEMENT	
	For Dominant Eye
Adjustment Steps	±0.25D

NEAR VISION ENHANCEMENT	
	For Non-Dominant Eye
Adjustment Steps	±0.25D

Plastic neutrality is established by purchasing credits from Plastic Bank. A credit represents the collection and conversion of one kilogram of plastic that may reach or be destined for waterways. CooperVision purchases credits equal to the weight of plastic in our one-day contact lens orders (effective from 1st January 2022). One-day contact lens plastic is determined by the weight of plastic in the blister, the lens and the secondary package (outer carton), including laminates, adhesives, and auxiliary inputs (e.g. ink). 1. CVI data on file 2020. Prospective, double-masked, bilateral, 1-week dispensing study with MyDay® daily disposable multifocal; n=104 habitual multifocal contact lens wearers. 2. CVI data on file 2021. Prospective, subject-masked, randomized, bilateral, two-week dispensing study at 5 US sites with MyDay® multifocal; n=58 habitual multifocal contact lens wearers. 3. CVI data on file 2021. Prospective, subject-masked, randomized, bilateral, two-week cross-over dispensing study at 5 US sites with MyDay® multifocal and DAILIES TOTAL1® Multifocal; n=58 habitual multifocal contact lens wearers. Packaging images for illustration purposes only.



MyDay[®] multifocal fitting guide

OptiExpert™ is available for mobile and tablet devices or as a web app, just scan the QR code of your choice

Benefits

- Simple to fit and an easy progression as the near requirement increases.¹
- Unsurpassed vision at all distances compared to other leading one-day multifocal contact lenses.^{2,3}
- Unsurpassed comfort compared to other leading one-day multifocal contact lenses.^{2,3}
- More prescription options than any other one-day multifocal contact lens.⁴
- UV-blocking.*
- Good handling.^{†§5}

Product specifications

Material	stenfilcon A
Water content	54%
Base curve	8.4mm
Diameter	14.2mm
Centre thickness (@ -3.00DS)	0.08mm
Dk/t (@ -3.00DS)	100
UV Blocker*	Yes
Power range	+8.00 to -10.00DS (0.25D steps) -10.50 to -12.00DS (0.50D steps)
ADD powers	LOW (+0.75 to +1.25D) MID (+1.50 to +1.75D) HIGH (+2.00 to +2.50D)
Multifocal design	Centre Near

Clinical tips




- Prescribe maximum plus power for binocular distance vision. **DO NOT OVER MINUS.**
- Use loose handheld lenses or flipper for over-refractions. **DO NOT USE A PHOROPTER.**
 - If distance vision needs to be enhanced, offer $\pm 0.25D$ to the dominant eye. If distance vision improves, check that near vision is maintained. Adjust the lens sphere power as applicable for the dominant eye. **DO NOT CHANGE ADD POWER.**
 - If near vision needs to be enhanced, offer $\pm 0.25D$ to the non-dominant eye. If near vision improves, check that distance vision is maintained. Adjust the lens sphere power as applicable for the non-dominant eye. **DO NOT CHANGE ADD POWER.**
- In some instances a refinement to the contralateral eye may improve distance/near vision.

CooperVision[®] Binocular Progressive System™

The Binocular Progressive System™ is designed to help presbyopes¹ achieve the vision experience they once had. It offers presbyopes¹ optimal visual acuity at all distances, even as their prescriptions and vision needs change.¹



Aquaform[®] Technology




Creates an optimised balance of high oxygen permeability and water content to provide high levels of breathability and moisture in a comfortable flexible lens.⁶

*Warning: UV-absorbing contact lenses are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, because they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed. † Vs. DAILIES TOTAL1[®] MULTIFOCAL and 1-DAY ACUVUE[®] MOIST MULTIFOCAL. ‡ Patients rated handling with MyDay[®] multifocal very good after a week of daily wear. § Mean rating 8.8/10. ¶ With myopia, hyperopia or emmetropia. 1. CVI data on file 2020. Prospective, double-masked, bilateral, 1-week dispensing study with MyDay[®] daily disposable multifocal; n=104 habitual multifocal contact lens wearers. 2. CVI data on file, 2020. Prospective, double-masked, randomized, bilateral, one-week cross-over dispensing study with MyDay[®] multifocal and 1 DAY ACUVUE[®] MOIST MULTIFOCAL; n=104 habitual MFCL wearers. 3. CVI data on file 2021. Prospective, subject-masked, randomized, bilateral, two-week cross-over dispensing study at 5 US sites with MyDay[®] multifocal and DAILIES TOTAL1[®] Multifocal; n=58 habitual MFCL wearers. 4. CooperVision data on file 2023. Based on prescription option combinations (sph + add) available across all daily disposable multifocal soft lenses from CVI, JIV, B+L and Alcon in UK, France, Germany and Italy Feb 2023. Cosmetic & Photochromatic CIs not included. Multiple base curve variants not included. 5. McParland M, Walsh K, Comoroda D, Sulley A. Wearer Experience and Eye Care Professional Acceptance with a 1 Day Multifocal Contact Lens. Clin Optom (Auckl). 2022 Nov 22;14:223-235. 6. CVI data on file, 2024.

clariti® 1 day multifocal fitting guide

Featuring smooth constant-powered progressive zones with four dedicated, discrete zones of stable power. Also featuring WetLoc® Technology

- Simple and successful fitting.¹
- 96% of patients successfully fitted with the first pairs of lenses.²
- 92% of patients agreed that clariti® 1 day multifocal met or exceeded their vision needs.³
- 98% of patients indicated that the comfort of the clariti® 1 day multifocal lenses met or exceeded expectation.⁴



Initial lens selection

- Step 1** – Using up-to-date spectacle prescription, determine ocular dominance using the **+1.00D blur** method.
- Step 2** – Determine spherical equivalent distance power (corrected for vertex distance).
- Step 3** – Select distance sphere power for each eye with ADD powers as indicated below.

FITTING GUIDE				
Spectacle Rx	Dominant Eye		Non-Dominant Eye	
ADD	Myope/ Emmetrope	Hyperope	Myope/ Emmetrope	Hyperope
+0.75 to +1.75D	BS LOW	BS LOW	BS LOW	BS +0.25D LOW
+2.00 to +2.25D	BS LOW	BS +0.25D LOW	BS +0.50D LOW	BS +0.50D LOW
>+2.25D	BS +0.25D LOW	BS +0.25D LOW	BS +0.25D HIGH	BS +0.25D HIGH

BS - Best sphere.

Plastic neutrality is established by purchasing credits from Plastic Bank. A credit represents the collection and conversion of one kilogram of plastic that may reach or be destined for waterways. CooperVision purchases credits equal to the weight of plastic in our one-day contact lens orders (effective from 1st January 2022). One-day contact lens plastic is determined by the weight of plastic in the blister, the lens and the secondary package (outer carton), including laminates, adhesives, and auxiliary inputs (e.g. ink). * Using OptiExpert™ tool. † At baseline approximately 5 mins after lens insertion. 1. CVI data on file 2013. Fitting assessment of c1dMF (n=22) and c1d (n=15) at 5 UK sites. Difference in average total chair time c1d 19 minutes and c1dMF 24 minutes. 2. Woods J et al. Validation of a multifocal contact lens online fitting app. BCLA poster presentation 2019. Retrospective analysis refraction data with OptiExpert™ (n=96 eyes) with Rx range +5.00 to -6.00DS, ≤ -1.00DC; ADDs +1.50 to +2.50D. 3. Based on retrospective analysis of 26 patients (52 eyes) with subjective refractions ranging between +5.00D to -6.00D and ≤ -1.00 D cylinder. 4. CVI data on file, 2019. Prospective, bilateral, subject-masked dispensing study with clariti® 1 day multifocal. N=48 habitual soft MFCL wearers at 2 sites in North America. Packaging images for illustration purposes only.

Vision assessment


- For best results, allow wearer to experience vision outside the testing room for 10-15minutes.
- Check vision with both eyes open and room lights on.
 - For **distance vision**, assess in surrounding environment under normal lighting conditions.
 - For **near vision**, assess using a mobile phone or other reading material.
- If vision acceptable, dispense trial lenses.
- If vision not acceptable, follow the lens optimisation steps described to the right.


Lens optimisation




Have patient keep both eyes open and optimise using handheld lenses or a flipper. **DO NOT USE A PHOROPTER.**

DO NOT CHANGE ADD POWER.



DISTANCE VISION ENHANCEMENT	
	For Dominant Eye
Adjustment Steps	±0.25D



NEAR VISION ENHANCEMENT	
	For Non-Dominant Eye
Adjustment Steps	±0.25D



clariti® 1 day multifocal fitting guide



OptiExpert™ is available for mobile and tablet devices or as a web app, just scan the QR code of your choice


Benefits

- Very good vision performance for near, intermediate or distance vision.¹
- Sustained high water content for excellent all-day comfort.
- 100% corneal oxygen consumption² for white, bright eyes.²
- UV-blocking.*
- Easy to handle.³

Product specifications


Material	somofilcon A
Water content	56%
Base curve	8.6mm
Diameter	14.1mm
Centre thickness (@ -3.00DS)	0.07mm
Dk/t (@ -3.00DS)	86
UV Blocker*	Yes
Power range	+5.00 to -6.00DS (0.25D steps)
ADD powers	LOW HIGH
Multifocal design	Centre Near

Clinical tips



- Prescribe maximum plus power for binocular distance vision. **DO NOT OVER MINUS.**
- Use loose handheld lenses or flipper for over-refractions. **DO NOT USE A PHOROPTER.**
 - If distance vision needs to be enhanced, offer $\pm 0.25D$ to the dominant eye. If distance vision improves, check that near vision is maintained. Adjust the lens sphere power as applicable for the dominant eye. **DO NOT CHANGE ADD POWER.**
 - If near vision needs to be enhanced, offer $\pm 0.25D$ to the non-dominant eye. If near vision improves, check that distance vision is maintained. Adjust the lens sphere power as applicable for the non-dominant eye. **DO NOT CHANGE ADD POWER.**
- In some instances a refinement to the contralateral eye may improve distance/near vision.

Dedicated discrete zones



Near Zone


Short Intermediate Zone

Long Intermediate Zone

Distance Zone

Unique dual Intermediate Vision Zones are designed to meet the lifestyle vision needs of many presbyopes.

WetLoc® Technology




WetLoc® Technology ensures the entire lens surface is continuously moist throughout your day. The built-in moisture is evenly spread throughout the lens and mimics the moisture dispersion of naturally healthy eyes.

Proclear® 1 day multifocal fitting guide

Designed with a single, aspheric profile for all add powers, and featuring PC technology™.

- Single power profile to preserve distance vision and limit visual compromise.
- Near boost in the non-dominant eye to optimise near and intermediate vision and maintain binocularity; accommodates patients up to +2.50 ADD.
- Easy adaptation for lens wearers through different stages of presbyopia.
- Material technology offers a natural resistance to dehydration.
- Soft contact lens material that is naturally biocompatible with the eye.



Initial lens selection

- Step 1** – Using up-to-date spectacle prescription, determine ocular dominance using the **+1.00D blur** method.
- Step 2** – Determine spherical equivalent distance power (corrected for vertex distance).
- Step 3** – Select distance sphere power for each eye with ADD powers as indicated below.

FITTING GUIDE		
Spectacle Rx	Dominant Eye	Non-Dominant Eye
ADD	Near Boost	Near Boost
Up to +1.00D	BS +0.50D	BS +0.50D
+1.25 to +2.50D	BS +0.50D	BS +1.25D

BS - Best sphere.

Plastic neutrality is established by purchasing credits from Plastic Bank. A credit represents the collection and conversion of one kilogram of plastic that may reach or be destined for waterways. CooperVision purchases credits equal to the weight of plastic in our one-day contact lens orders (effective from 1st January 2022). One-day contact lens plastic is determined by the weight of plastic in the blister, the lens and the secondary package (outer carton), including laminates, adhesives, and auxiliary inputs (e.g. ink). Packaging images for illustration purposes only.

Vision assessment


- For best results, allow wearer to experience vision outside the testing room for 10-15 minutes.
- Check vision with both eyes open and room lights on.
 - For **distance vision**, assess in surrounding environment under normal lighting conditions.
 - For **near vision**, assess using a mobile phone or other reading material.
- If vision acceptable, dispense trial lenses.
- If vision not acceptable, follow the lens optimisation steps described to the right.


Lens optimisation




Have patient keep both eyes open and optimise using handheld lenses or a flipper.
DO NOT USE A PHOROPTER.

DO NOT CHANGE ADD POWER.



DISTANCE VISION ENHANCEMENT	
	For Dominant Eye
Adjustment Steps	±0.25D



NEAR VISION ENHANCEMENT	
	For Non-Dominant Eye
Adjustment Steps	±0.25D



Proclear® 1 day multifocal fitting guide



OptiExpert™ is available for mobile and tablet devices or as a web app, just scan the QR code of your choice

Benefits

- Excellent distance and intermediate vision and very good near vision.*
- Retains up to 96% of their original water content, even after 12 hours of wear.¹
- May help address eye dryness when wearing contact lenses.

Product specifications



Material	omafilcon A
Water content	60%
Base curve	8.7mm
Diameter	14.2mm
Centre thickness (@ -3.00DS)	0.09mm
Dk/t (@ -3.00DS)	28
UV Blocker	No
Power range	-0.25 to -6.00DS (0.25D steps) -6.50 to -10.00DS (0.50D steps) +0.25 to +6.00DS (0.25D steps)
ADD powers	Single power profile +1.50D (Could fit adds up to +2.50)
Multifocal design	Centre Near

Clinical tips



- Prescribe maximum plus power for binocular distance vision. **DO NOT OVER MINUS.**
- Use loose handheld lenses or flipper for over-refractions. **DO NOT USE A PHOROPTER.**
 - If distance vision needs to be enhanced, offer $\pm 0.25D$ to the dominant eye. If distance vision improves, check that near vision is maintained. Adjust the lens sphere power as applicable for the dominant eye. **DO NOT CHANGE ADD POWER.**
 - If near vision needs to be enhanced, offer $\pm 0.25D$ to the non-dominant eye. If near vision improves, check that distance vision is maintained. Adjust the lens sphere power as applicable for the non-dominant eye. **DO NOT CHANGE ADD POWER.**
- In some instances a refinement to the contralateral eye may improve distance/near vision.

PC Technology™



CooperVision's exclusive PC Technology™ creates an ideal lens for presbyopic patients and may help those with dry eyes.

Based on molecular chemistry, PC Technology™ creates a unique lens material in which the phosphorylcholine (PC) molecules attract and bind water to the surface, creating a shield that keeps the lenses clean and functioning properly.

The PC molecules also help the lenses remain hydrated, which in turn, help them feel moist and comfortable all day long.

* in high contrast situations.

1. Hall B, Jones S, Young G, Coleman S. The on-eye dehydration of proclear compatibles lenses. CLAO J. 1999 Oct;25(4):233-7. PMID: 10555740.

Biofinity® multifocal fitting guide




Featuring **Balanced Progressive® Technology** and **Aquaform® Technology**

- Two different optical designs (D and N) to enhance and provide exceptional vision at all distances – near, intermediate and far.
- Lens design is further optimised for each sphere and ADD power.
- High oxygen permeability and naturally wettable.
- 98% of patients successfully fitted using two pairs of lenses or fewer.^{†1}



Initial lens selection

- Step 1** – Using up-to-date spectacle prescription, determine ocular dominance using the **+1.00D blur** method.
- Step 2** – Determine spherical equivalent distance power (corrected for vertex distance).
- Step 3** – Select distance sphere power for each eye with ADD powers as indicated below.

FITTING GUIDE		
 Spectacle Rx	 Dominant Eye	 Non-Dominant Eye
ADD*	DESIGN	DESIGN
+1.00D	D	D
+1.50D	D	D
+2.00D	D	N
+2.50D	D	N

D refers to a centre Distance design. N refers to a centre Near design. *Always round down to the nearest available ADD. † with the OptiExpert™ trial lens. 1. CVI data on file 2019. Retrospective analysis; N=55 subjects (110 eyes); DV Rx +1.25D to -3.25D, ADD powers +1.25 to +2.50DS. Packaging images for illustration purposes only.

Vision assessment


- For best results, allow wearer to experience vision outside the testing room for 10-15 minutes.
- Check vision with both eyes open and room lights on.
 - For **distance vision**, assess in surrounding environment under normal lighting conditions.
 - For **near vision**, assess using a mobile phone or other reading material.
- If vision acceptable, dispense trial lenses.
- If vision not acceptable, follow the lens optimisation steps described to the right.


Lens optimisation




Have patient keep both eyes open and optimise using handheld lenses or a flipper. **DO NOT USE A PHOROPTER.**

DO NOT CHANGE ADD POWER.



DISTANCE VISION ENHANCEMENT	
	For Dominant Eye
Adjustment Steps	±0.25D



NEAR VISION ENHANCEMENT	
	For Non-Dominant Eye
Adjustment Steps	±0.25D



Biofinity® multifocal fitting guide



OptiExpert™ is available for mobile and tablet devices or as a web app, just scan the QR code of your choice

Benefits

- Allows for both a simplified and flexible fitting for presbyopic patients.
- Allows for an individualised fitting for each wearer and eye.
- Superior visual clarity at all distances.
- Incredible all-day comfort.
- Deliver a high level of oxygen to your presbyopic patients' eyes.¹

Product specifications


Material	comfilcon A
Water content	48%
Base curve	8.6mm
Diameter	14.0mm
Centre thickness (@ -3.00DS)	0.09mm
Dk/t (@ -3.00DS)	142
UV Blocker	No
Power range	+6.00 to -6.00DS (0.25D steps) -6.50 to -10.00DS (0.50D steps)
ADD powers	+1.00, +1.50, +2.00, +2.50D D&N
Multifocal design	Centre Distance and Centre Near

Clinical tips




- Prescribe maximum plus power for binocular distance vision. **DO NOT OVER MINUS.**
- Use loose handheld lenses or flipper for over-refractions. **DO NOT USE A PHOROPTER.**
 - If distance vision needs to be enhanced, offer $\pm 0.25D$ to the dominant eye. If distance vision improves, check that near vision is maintained. Adjust the lens sphere power as applicable for the dominant eye. **DO NOT CHANGE ADD POWER.**
 - If near vision needs to be enhanced, offer $\pm 0.25D$ to the non-dominant eye. If near vision improves, check that distance vision is maintained. Adjust the lens sphere power as applicable for the non-dominant eye. **DO NOT CHANGE ADD POWER.**
- In some instances a refinement to the contralateral eye may improve distance/near vision.

Balanced Progressive® Technology



Two different optical designs utilise the processing power of the visual cortex to enhance vision. Optimised for exceptional vision at all distances – near, intermediate and far. Lens design is further optimised for each sphere and ADD power.

Aquaform® Technology




Creates an optimised balance of high oxygen permeability and water content to provide high levels of breathability and moisture in a comfortable flexible lens.¹

Biofinity® toric multifocal fitting guide

Featuring Optimised Toric Lens Geometry™, Balanced Progressive® Technology and Aquaform® Technology




- A stable fit with remarkable visual performance.
- 93% of patients successfully fitted with the first pairs of lenses.¹
- Optimised Toric Lens Geometry™ provides excellent stability, predictable orientation, and consistent visual acuity.
- Balanced Progressive® Technology optimised for each sphere and ADD power.
- Patients rated Biofinity® toric multifocal lenses 8 out of 10 for comfort during the day.²



Initial lens selection

- Step 1** – Using up-to-date spectacle prescription, determine ocular dominance using the **+1.00D blur** method.
- Step 2** – Determine the sphere and cylinder power and axis rounding to the nearest 5° if necessary.
- Step 3** – **(optional)** Use Biofinity® toric fit set to confirm the toric trial lens parameter. Adjust axis based on rotation, rounding to the nearest 5° if necessary.

Select distance sphere power for each eye with ADD powers as indicated below:

FITTING GUIDE		
 Spectacle Rx	 Dominant Eye	 Non-Dominant Eye
ADD*	DESIGN	DESIGN
+1.00D	D	D
+1.50D	D	D
+2.00D	D	N
+2.50D	D	N

Vision assessment


- For best results, allow wearer to experience vision outside the testing room for 10–15 minutes.
- Check vision with both eyes open and room lights on.
 - For **distance vision**, assess in surrounding environment under normal lighting conditions.
 - For **near vision**, assess using a mobile phone or other reading material.
- Assess toric orientation and adjust axis if required before refining multifocal powers.
- If vision acceptable, dispense trial lenses.
- If vision not acceptable, follow the lens optimisation steps described to the right.


Lens optimisation




Have patient keep both eyes open and optimise using handheld lenses or a flipper.
DO NOT USE A PHOROPTER.

DO NOT CHANGE ADD POWER.

DISTANCE VISION ENHANCEMENT	
	For Dominant Eye
Adjustment Steps	±0.25D

NEAR VISION ENHANCEMENT	
	For Non-Dominant Eye
Adjustment Steps	±0.25D

Biofinity® toric multifocal fitting guide



Benefits

- Provides predictable, consistent fitting for astigmatic presbyopic patients.
- Allows for an individualised fitting for each wearer and eye.
- Offers the excellent vision performance your astigmatic presbyopic patients expect and deserve.
- Allows for clear vision – near, far and in-between.
- Stays comfortable all day.
- Delivers plenty of oxygen to your patients' eyes.¹

Product specifications

Material	comfilcon A
Water content	48%
Base curve	8.7mm
Diameter	14.5mm
Centre thickness (@ -3.00DS)	0.11mm
Dk/t (@ -3.00DS)	116
UV Blocker	No
Power range	-10.00 to -6.50DS (0.50D steps) -6.00 to +6.00DS (0.25D steps) +6.50 to +10.00DS (0.50D steps)
Cylinder powers	-0.75 to -5.75 (0.50 steps)
ADD powers	+1.00D, +1.50D, +2.00D, +2.50D
Axes	5° to 180° (in 5° steps)
Multifocal design	Centre Distance and Centre Near


Made to order.
1. CVI Data on file 2024.

Clinical tips




- Prescribe maximum plus power for binocular distance vision. **DO NOT OVER MINUS.**
- Use loose handheld lenses or flipper for over-refractions. **DO NOT USE A PHOROPTER.**
 - If distance vision needs to be enhanced, offer $\pm 0.25D$ to the dominant eye. If distance vision improves, check that near vision is maintained. Adjust the lens sphere power as applicable for the dominant eye. **DO NOT CHANGE ADD POWER.**
 - If near vision needs to be enhanced, offer $\pm 0.25D$ to the non-dominant eye. If near vision improves, check that distance vision is maintained. Adjust the lens sphere power as applicable for the non-dominant eye. **DO NOT CHANGE ADD POWER.**
- In some instances a refinement to the contralateral eye may improve distance/near vision.

Optimised Toric Lens Geometry™




The smooth, continuous ballast maximises comfort during lid interaction and a stable fit, whilst the large toric optic zone helps provide excellent vision performance. A unique selection of back surface curvatures optimise on-eye movement to support stable lens positioning, for clear vision a predictable fit, and comfort.

Balanced Progressive® Technology



Two different optical designs utilise the processing power of the visual cortex to enhance vision. Optimised for exceptional vision at all distances – near, intermediate and far. Lens design is further optimised for each sphere and ADD power.

Aquaform® Technology



Creates an optimised balance of high oxygen permeability and water content to provide high levels of breathability and moisture in a comfortable flexible lens.¹

Proclear® multifocal/multifocal XR fitting guide




Featuring Balanced Progressive® Technology and PC Technology™

- Multiple zones of vision correction.
- Natural resistance to dehydration.
- Biocompatible with the eye.



Initial lens selection

- Step 1** – Using up-to-date spectacle prescription, determine ocular dominance using the **+1.00D blur** method.
- Step 2** – Determine spherical equivalent distance power (corrected for vertex distance).
- Step 3** – Select distance sphere power for each eye with ADD powers as indicated below.

FITTING GUIDE		
 Spectacle Rx	 Dominant Eye	 Non-Dominant Eye
ADD*	DESIGN	DESIGN
+1.00D	D	D
+1.50D	D	D
+2.00D	D	N
+2.50D	D	N

D refers to a centre Distance design. N refers to a centre Near design. *Always round down to the nearest available ADD.
Packaging images for illustration purposes only.

Vision assessment


- For best results, allow wearer to experience vision outside the testing room for 10-15 minutes.
- Check vision with both eyes open and room lights on.
 - For **distance vision**, assess in surrounding environment under normal lighting conditions.
 - For **near vision**, assess using a mobile phone or other reading material.
- If vision acceptable, dispense trial lenses.
- If vision not acceptable, follow the lens optimisation steps described to the right.


Lens optimisation




Have patient keep both eyes open and optimise using handheld lenses or a flipper.
DO NOT USE A PHOROPTER.

DO NOT CHANGE ADD POWER.

DISTANCE VISION ENHANCEMENT	
	For Dominant Eye
Adjustment Steps	±0.25D

NEAR VISION ENHANCEMENT	
	For Non-Dominant Eye
Adjustment Steps	±0.25D

Proclear® multifocal/multifocal XR fitting guide



Benefits

- Allows for clear vision – near, far, and in-between.
- All-day comfort.¹
- Can stay 96% hydrated throughout the day, even after 12 hours of wear.
- Attracts and maintains moisture, reducing dryness.
- Extensive parameter range for presbyopic patients.

Product specifications


Material	omafilcon B
Water content	62%
Base curve	8.7mm
Diameter	14.4mm
Centre thickness (@ -3.00DS)	0.14mm
Dk/t (@ -3.00DS)	17
UV Blocker	No
Power range	+6.50 to -6.50DS (0.25D steps) +7.00 to +20.00DS (0.50D steps) -7.00 to -20.00DS (0.50D steps)
ADD powers	+1.00 to +2.50DS (0.50D steps) +3.00 to +4.00DS (0.50D steps)
Multifocal design	Centre Distance and Centre Near

Clinical tips



- Prescribe maximum plus power for binocular distance vision. **DO NOT OVER MINUS.**
- Use loose handheld lenses or flipper for over-refractions. **DO NOT USE A PHOROPTER.**
 - If distance vision needs to be enhanced, offer $\pm 0.25D$ to the dominant eye. If distance vision improves, check that near vision is maintained. Adjust the lens sphere power as applicable for the dominant eye. **DO NOT CHANGE ADD POWER.**
 - If near vision needs to be enhanced, offer $\pm 0.25D$ to the non-dominant eye. If near vision improves, check that distance vision is maintained. Adjust the lens sphere power as applicable for the non-dominant eye. **DO NOT CHANGE ADD POWER.**
- In some instances a refinement to the contralateral eye may improve distance/near vision.

PC Technology™



Creates a lens material that contains phosphorylcholine (PC) molecules, which bind with water molecules in natural tears to create a 'shield' of water around the lens.

Proclear® multifocal toric fitting guide




Featuring Balanced Progressive® Technology and PC Technology™

- Multiple zones of vision correction.
- Natural resistance to dehydration.
- Biocompatible with the eye.



Initial lens selection

- Step 1** – Using up-to-date spectacle prescription, determine ocular dominance using the **+1.00D blur** method.
- Step 2** – Determine spherical equivalent distance power and axis, rounding to the nearest 5° if necessary.
- Step 3** – **(optional)** Use Proclear® toric fit set to confirm the toric trial lens axes. Adjust axis based on rotation, rounding to the nearest 5° if necessary.

FITTING GUIDE		
 Spectacle Rx	 Dominant Eye	 Non-Dominant Eye
ADD*	DESIGN	DESIGN
+1.00D	D	D
+1.50D	D	D
+2.00D	D	N
+2.50D	D	N

D refers to a centre Distance design. N refers to a centre Near design. *Always round down to the nearest available ADD.
Packaging images for illustration purposes only.

Vision assessment


- For best results, allow wearer to experience vision outside the testing room for 10–15 minutes.
- Check vision with both eyes open and room lights on.
 - For **distance vision**, assess in surrounding environment under normal lighting conditions.
 - For **near vision**, assess using a mobile phone or other reading material.
- Assess toric orientation and adjust axis if required before optimising multifocal powers.
- If vision acceptable, dispense trial lenses.
- If vision not acceptable, follow the lens optimisation steps described to the right.


Lens optimisation




Have patient keep both eyes open and optimise using handheld lenses or a flipper.
DO NOT USE A PHOROPTER.

DO NOT CHANGE ADD POWER.



DISTANCE VISION ENHANCEMENT	
	For Dominant Eye
Adjustment Steps	±0.25D



NEAR VISION ENHANCEMENT	
	For Non-Dominant Eye
Adjustment Steps	±0.25D



Proclear® multifocal toric fitting guide



OptiExpert™ is available for mobile and tablet devices or as a web app, just scan the QR code of your choice

Benefits

- Vision correction for astigmatic presbyopes, including higher prescriptions.
- Allows for clear vision – near, far, and in-between.
- All-day comfort.¹
- Can stay 96% hydrated throughout the day, even after 12 hours of wear.
- Attracts and maintains moisture, reducing dryness.

Product specifications

Material	omafilcon B
Water content	62%
Base curve	8.4, 8.8mm
Diameter	14.4mm
Centre thickness (@ -3.00DS)	Varies by Rx
Dk/t (@ -3.00DS)	Varies by Rx
UV Blocker	No
Power range	+6.50 to -6.50DS (0.25D steps) +7.00 to +20.00DS (0.50D steps) -7.00 to -20.00DS (0.50D steps)
Cylinder powers	-0.75 to -5.75DC (0.50D steps)
ADD powers	+1.00 to +4.00DS (0.50D steps)
Axes	5° to 180° (in 5° steps)
Multifocal design	Centre Distance and Centre Near

Made to order.


1. Hall B et al. The On-Eye Dehydration of Proclear Compatibles Lenses. CLAO (1999); 25(4): 233-237.

Clinical tips



- Prescribe maximum plus power for binocular distance vision. **DO NOT OVER MINUS.**
- Assess toric orientation and adjust axis if required before refining multifocal powers.
- Use loose handheld lenses or flipper for over-refractions. **DO NOT USE A PHOROPTER.**
 - If distance vision needs to be enhanced, offer $\pm 0.25D$ to the dominant eye. If distance vision improves, check that near vision is maintained. Adjust the lens sphere power as applicable for the dominant eye. **DO NOT CHANGE ADD POWER.**
 - If near vision needs to be enhanced, offer $\pm 0.25D$ to the non-dominant eye. If near vision improves, check that distance vision is maintained. Adjust the lens sphere power as applicable for the non-dominant eye. **DO NOT CHANGE ADD POWER**
- In some instances a refinement to the contralateral eye may improve distance/near vision.


PC Technology™



Creates a lens material that contains phosphorylcholine (PC) molecules, which bind with water molecules in natural tears to create a 'shield' of water around the lens.

OptiExpert™ v2.0 app.

Smarter than ever, so you can do more.



Instantly converts (virtually) any spectacle prescription into a contact lens prescription.

OptiExpert™ v2.0 helps make contact lens selection even easier.*

Smart Prescription Calculator: Instantly converts virtually any sphere, toric, multifocal or toric multifocal spectacle prescription into a contact lens prescription.

Efron Grading Scales: Readily accessible clinical reference.*

Oxygen Profiles: Effective clinical tool to simplify patient communications.*



OptiExpert™ is available for mobile and tablet devices or as a web app.



*OptiExpert™ is an educational, reference and information tool for Eye Care Professionals. Eye Care Professionals may choose to use the app in connection with their own patient evaluation but it is not intended to be relied upon for clinical decision-making. OptiExpert™ is not intended as and does not constitute medical or optometric advice nor is it intended to replace the patient evaluation performed by an Eye Care Professional.



CooperVision®

www.coopervision.co.uk

