

MyDay® daily disposable multifocal fitting guide.

Featuring the CooperVision® Binocular Progressive System.™



- Designed for simplicity of fit.
- Highly effective and successful^{1,2} fitting approach.
- Supported by OptiExpert.™

Initial lens selection.

Step 1 Using up-to-date spectacle prescription, determine spherical equivalent distance power (corrected for vertex distance).

Step 2 Determine eye dominance with **+1.00D** blur method.

Step 3 Select distance sphere power for each eye with add powers as indicated below.

INITIAL CONTACT LENS SELECTION		
Spectacle Rx Add	Dominant Eye	Non-Dominant Eye
+0.75D to +1.25D	LOW	LOW
+1.50D to +1.75D	LOW	MED
+2.00D to +2.50D	LOW	HIGH

Vision assessment.

- For best results, allow wearer to experience vision outside the testing room for 10-15 minutes.
- Check vision with both eyes open and room lights on.
 - For **distance vision**, assess in surrounding environment under normal lighting conditions.
 - For **near vision**, assess using a mobile phone or other reading material.
- If vision acceptable, dispense trial lenses.
- If vision not acceptable, follow the lens optimisation steps described to the right.

Lens optimisation.



Have patient keep both eyes open and optimise using handheld lenses or a flipper; do not use a phoropter.

DO NOT CHANGE ADD POWER.

DISTANCE VISION ENHANCEMENT	
	For Dominant Eye
Adjustment Steps	±0.25D
NEAR VISION ENHANCEMENT	
	For Non-Dominant Eye
Adjustment Steps	±0.25D

1. CVI data on file 2020. Prospective, double-masked, bilateral, one-week dispensing study UK with MyDay® daily disposable multifocal; n=104 habitual multifocal contact lens wearers;

2. CVI data on file 2021. Prospective, subject-masked, randomised, bilateral, two-week dispensing study at five US sites with MyDay® daily disposable multifocal; n=58 habitual multifocal contact lens wearers.
MyDay MF Fitting Guide

MyDay® daily disposable multifocal fitting guide.



OptiExpert™ is available for mobile and tablet devices or as a web app.



Product specifications.

MyDay® daily disposable multifocal.

Base Curve	8.4 mm
Diameter	14.2 mm
Power Range	+8.00D to -10.00D (0.25D steps) -10.50D to -12.00D (0.50D steps)
Add Power	Low (+0.75D to +1.25D spectacle Rx add) Med (+1.50D to +1.75D spectacle Rx add) High (+2.00D to +2.50D spectacle Rx add)
Material	stenfilcon A
Dk/t (at -3.00D)	100
Centre thickness	0.08mm
UV Blocking*	85% UVA/96% UVB

*Warning: UV-absorbing contact lenses are not substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, because they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed.

1. CVI data on file, 2020. Prospective, double-masked, bilateral, one-week dispensing study with MyDay® daily disposable multifocal; n=104 habitual multifocal contact lens wearers.

MyDay MF Fitting Guide

©2022 CooperVision 02/22.



Clinical tips.



- Prescribe maximum plus power for binocular distance vision; do not over minus.
- Use loose handheld lenses or flipper for over-refractions. Do not use a phoropter.
 - If distance vision needs to be enhanced, offer $\pm 0.25D$ to the dominant eye. If distance vision improves, check that near vision is maintained. Adjust the lens sphere power as applicable for the dominant eye. **DO NOT CHANGE ADD POWER.**
 - If near vision needs to be enhanced, offer $\pm 0.25D$ to the non-dominant eye. If near vision improves, check that distance vision is maintained. Adjust the lens sphere power as applicable for the non-dominant eye. **DO NOT CHANGE ADD POWER.**

CooperVision® Binocular Progressive System.™

The Binocular Progressive System™ is designed to help presbyopes achieve the vision experience they once had. It offers presbyopes optimal visual acuity at all distances,¹ even as their prescriptions and vision needs change.



For illustrative purposes only.