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

Australia

SYMBOLS KEY:

The following symbols may appear on the label or carton.

SYMBOL	DEFINITION
	Caution: By law, this device must be prescribed by a registered practitioner
	Caution / See Instructions for Wearers
	Use by Date (expiration date)
	Batch Code
	Sterile using Steam Heat
	UV Blocking
	Manufacturer
	Authorized representative in the European Community
	Medical device
	Do not use if package is damaged
	Consult instructions for use
	Do not re-use

Soft contact lenses in buffered saline.

Always wash and dry your hands before handling your lenses.

Do not use if packaging is not sealed.

Contact lenses are not to be shared.

Lens Insertion

1. Before insertion, ensure the lens is not inside out. Place on your index finger to verify the shape (Fig. 1 & 2).
2. Using the middle finger, hold your lower lid down. With the index finger of the other hand, raise upper lid. Place the lens onto the eye (Fig. 3). Repeat for the other eye.

Lens Removal

Separate the lids with the middle finger of each hand. Look up and with the index finger of the lower hand; slide the lens down onto the white of the eye. Pinch the lens and lift out with the thumb and index finger (Fig. 4). Repeat for the other eye.

Disinfection/Soaking

For lenses worn more than once clean and disinfect them following strictly the instruction on your solution bottle or as directed by your eye care practitioner.

N.B. Daily disposable lenses are intended to be worn only once and then discarded. If your lenses are intended for daily disposable wear, DO NOT wear them more than once. Daily disposable lenses are not designed to be cleaned or disinfected and worn more than one time. Daily disposable lenses worn more than one time may become damaged and could put you at great risk of contact lens care and handling related eye infections or other problems.

If the lenses dry out:

Re-wet with disinfecting, soaking or rinsing solution then follow the Disinfection/Soaking step above

Indications

For the correction of ametropia.

Precautions:

1. For all types of lenses, your eye care practitioner will recommend a lens wearing, cleaning and replacement regime. Do not vary these. If in doubt seek the advice of your eye care practitioner.
2. If wearing lenses on an extended wear basis, remove and dispose of your lenses the evening before you are due to replace them, inserting new lenses the following day.
3. For lenses worn more than once, always clean and disinfect your lenses as instructed, after lens removal.
4. If your eyes become red or irritated or you experience any pain or unexpected change in vision, remove your lenses immediately and consult your eye care practitioner.
5. Avoid contaminating lenses with make-up, aerosol sprays or wearing them near noxious fumes.
6. Replace damaged, contaminated or scratched lenses immediately.
7. Avoid long nails, which harbour bacteria and could damage the eye or lens. A scratched eye is vulnerable to infection.
8. Keep your lenses out of the reach of children. *
9. Only use your lenses for water sports if watertight goggles are worn.
10. Contact with unsterile water (e.g. from the tap, bottled, bath, pool, or water sports) increases the risk of microbial infection.

* exception MiSight® 1 day contact lenses

Warnings:

1. **DO NOT** sleep in your lenses unless your eye care practitioner has advised it is safe to do so.
2. **DO NOT** wear your lenses beyond the period recommended by your eye care practitioner.
3. **DO NOT** use household products (e.g. disinfectants) on your lenses.
4. **DO NOT** wear your lenses when using eye drops (unless advised to do so by a medical practitioner).
5. **DO NOT** use saliva to wet your lenses.
6. **DO NOT** rinse your lenses or lens case with tap water.

Contraindications (reasons not to use)

Do not use contact lenses when any of the following conditions exist.

- Acute and subacute inflammation or infection of the anterior chamber of the eye.
- Any eye disease, injury, or abnormality that affects the cornea, conjunctiva, or eyelids.
- Corneal hypoesthesia (reduced corneal sensitivity), if not aphakic.
- Severe insufficiency of lacrimal secretion (dry eyes).
- Any systemic disease that may affect the eye or be exaggerated by wearing contact lenses.
- Allergic reactions of ocular surfaces or adnexa that may be induced or exaggerated by wearing contact lenses or use of contact lens solutions.
- Any active corneal infection (bacterial, fungal, protozoal or viral).
- If eyes become red or irritated.
- Use of a medication that is contraindicated, including eye medications.

Sponsor:

CooperVision Australia Pty Ltd
115, Hayward Avenue, Torrensville
SA 5031, Australia

Legal Manufacturer:

CooperVision Manufacturing Limited,
South Point, Hamble, Southampton,
SO31 4RF, UK

  CooperVision
Hamble, SO31 4RF, UK
0123   www.coopervision.com

 MD

Medical
Device

 EC REP

CooperVision CL Kft.
Gorcsev Iván street 7. Building C.
2360 Gyál, Hungary

STERILE 