


clariti® 1 day toric

High breathability¹, comfort and great value[†] in convenient daily disposable contact lenses.



CooperVision®
Live Brightly.®



If you need contact lenses for astigmatism, are health-conscious and are looking for affordable[†], high quality comfortable soft lenses that will keep you feeling good about your budget and your eyes, look no further than clariti® 1 day daily toric lenses—from CooperVision®, a world leader in toric contact lenses.

With clariti® 1 day toric, you can be confident that your eye health is supported because you replace each lens daily, it has a built-in UV blocker* and it's made with a healthy lens material that delivers 100% of the oxygen your eyes need.^{††1} While wearing clariti® 1 day toric contact lenses, you'll enjoy great vision², dependable stability and a comfortable lens-wearing experience, all at a price that'll make you smile.[†]

Replacement Schedule



Daily

Corrects



Astigmatism

Features & Benefits

- CooperVision® toric design provides a consistently stable fit you can depend on.
- Soft silicone hydrogel lens material with high breathability helps keep your eyes healthy and white.^{†1}
- UV blocker helps protect your eyes against the transmission of ultraviolet rays, helping support your eye health.*
- WetLoc® Technology and high water content help keep your contact lenses moist, supporting excellent all-day wearer comfort.
- clariti® 1 day has all the benefits of a silicone hydrogel daily disposable without the price premium.[†]

* Warning: UV-absorbing contact lenses are not substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, because they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed.

† During daily wear.

†† With higher oxygen permeability than hydrogel materials, silicone hydrogel contact lenses minimise or eliminate hypoxia-related signs and symptoms during lens wear.

† Based on CooperVision's recommended retail price. Retail price is solely at the discretion of the retailer at all times.

The above is not meant to influence the price any retailer charges.

References:

1. Brennan NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygenation During Contact Lens Wear. *Optom Vis Sci* 2005. 2005;82(6):467-472.
2. CVI data on file, 2018. Clinical performance study with clariti 1 day toric in UK, n=37. High contrast VA.

©2021 CooperVision

clariti®, CooperVision®, and WetLoc® Technology are registered trademarks of The Cooper Companies, Inc., and its subsidiaries.